



AED Maintenance Checklist

AED Location	
Make/ Model	
Serial Number	

	Week	1	2	3	4	5	6	7	8	9	10
AED Location											
Visible & Unobstructed	✓										
AED Sign Displayed	✓										
AED Condition											
Battery Status/ Date	✓										
Indicator Light OK	✓										
AED Exterior OK	✓										
AED Components											
Adult Pads	✓										
Child Pads or Switch *	✓										
Check expiry date on pads	✓										
CPR Mask	✓										
Clothing Shears	✓										
Safety Razor	✓										
Gloves	✓										
Initial											
Date											



To Order Supplies please visit www.hearts.ie/shop